

## Specifications Multicheck

### Function On Control Acid

Meas	20g	300mg/dL
(1.1)	~	33.3 mmol/L
100	~	400 mg/dL
(2.5)	~	10.4 mmol/L
3	~	20 mg/dL
(179)	~	1190 $\mu$ mol/L

Calibration Plasma equivalent  
 Test time 6 seconds 5 seconds 6 seconds

Memory capacity 200 test results 50 test results 10 test results

Operating condition  
 Temperature, humidity 57.2 ~ 104°F (14 ~ 40°C)  
 ≤ 85% Relative Humidity  
 Meter storage condition (*Transportation condition*)  
 Temperature, humidity 14 ~ 140°F (-10 ~ 60°C)  
 ≤ 95% Relative Humidity

Sample volume ≤  
 45  $\mu$ L

Sample type Finger capillary whole blood  
 Hematocrits range 30 ~ 50%  
 Power supply 1.5V (AAA) x 2  
 Meter Dimension HxWxD (mm) 88 x 64 x 22  
 Display LCD display (35 x 45 mm)  
 Weight 59 grams, without batteries  
 Life of battery More than 1.000 times  
 Technology Used Electrode-based Biosensor

- \* There will be 100 test results memories for dual-function meter or mono-function meter.
- \* There will be 50 test results memories for triple-function meter.
- \* The meter has been safety and EMC tested and approved according to the requirements for EN 61010-1/EN 61010-2-101/EN60601-1-2/EN61326

## Contents



Each pack contains:

- 1 Kernel Multicheck Monitoring System
- Kernel MultiCheck user's manual in Greek or English
- Quick user manual
- 1 Puncturer (Lancing device)
- Log book
- Pouch
- Lancets
- Check strip
- 2 AAA batteries
- Test strips (optional)
- Warranty card

### Instructions

The advantage of **Multicheck Monitoring System by Kernel** is that it is very easy to use and it allows a very short time to learn the exact levels of **glucose**

,  
**cholesterol**  
or  
**uric acid**  
in your blood.



#### Step 1

Take a blood sample from your finger.



#### Step 2

Place one drop of blood on the test strip and insert it in the Meter.



### Step 3

Results for Glucose: 6 seconds

Results for Cholesterol: 150 seconds

Results for Uric Acid: 6 seconds. Get your measurement set you and make your life easier!

## Customer Service

Our customer service includes:

(1) Our customers are entitled to **free repair and replacement parts**. If you have any problems with this product, please contact your local distributor / dealer for more information. Excluded are the damages caused by misuse or accidents (such as breaking, moisture, unauthorized repair, disassembly and natural disasters).

(2) The company Kernel Int'l Corp. assumes no responsibility for any accidents and exacerbation of disease resulting from the incorrect use of the device by the buyer / user without the guidance of professionals (eg medical staff). Consumers can claim compensation in this case.

(3) Please read the instruction manual thoroughly **before** using the meter Multicheck. If you have any questions, please contact us on (+30) 2510 391 910.